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A small, stylized handwritten signature in black ink.

An official circular stamp with illegible text around the perimeter. Overlaid on the stamp is handwritten text in Cyrillic: "Законно" (Legal), "Ма. 06.2017" (May 06, 2017), and "РА-2017-110".

A small handwritten mark or signature consisting of a circle with a horizontal line extending to the right.

A small handwritten mark consisting of a vertical line with a horizontal tick at the top.

A small handwritten mark consisting of a circle with a vertical line extending downwards.

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# CERTIFICATE

No. Q1N 14 07 68207 025

**Holder of Certificate:** **MP Biomedicals Asia Pacific Pte Ltd**  
 2 Pioneer Place  
 Singapore 627885  
 SINGAPORE

**Facility(ies):** MP Biomedicals Asia Pacific Pte Ltd  
 2 Pioneer Place, Singapore 627885, SINGAPORE

**Certification Mark:**

*Handwritten signatures and stamps:*  
 J. Wirth  
 Director  
 TUV SUD

**Scope of Certificate:** **Design, Production, Storage and Distribution of In-vitro Diagnostics for the detection of infectious diseases and pregnancy; Distribution of Western Blot and ELISA Instruments; Production and Purification of Antigen and Protein**

**Applied Standard(s):**

EN ISO 13485:2012 + AC:2012  
 Medical devices - Quality management systems - Requirements for regulatory purposes (ISO 13485:2003 + Cor. 1:2009)  
 DIN EN ISO 13485:2012

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). See also notes overleaf.

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
Diagnostics Division  
29525 Fountain Parkway  
Solon, OH 44139  
<http://www.mpbio.com>

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Customer support and technical information can be obtained at local MP Biomedicals offices.

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F. Extracting Reagent, **REAG EXT** Catalog No. 06B257176, contains 1.0N sodium hydroxide with organic extracting enhancer and yellow dye. Volume: >10 mL/bottle. One bottle /100 tube kit, 2 bottles/200 tube kit. **Storage:** Refrigerate at 2-8°C. **Stability:** 24 months from date of manufacture.

**CAUTION:** Avoid contact with eyes, skin and clothing. Take appropriate safety precautions when pipetting this reagent.

**\*WARNING: HANDLE AS IF CAPABLE OF TRANSMITTING INFECTION:** Source material from which this product was derived was found nonreactive for HBsAg and negative for HIV antibody when tested with licensed reagents. No known test method can offer assurance that products derived from human blood will not be infectious. Refer to CDC/NIH Biosafety in Microbiological and Biomedical Laboratories Publication (IHS Publication No. [CDC] 84-8395).

**WARNING: CONTAINS RADIOACTIVE MATERIAL**

This MP Biomedicals Radioassay Kit contains <3 microcuries (111 kilobecquerels) of [<sup>125</sup>I] and <1.5 microcuries (55.5 kilobecquerels) of [<sup>57</sup>Co] per vial of tracer. This radioactive material may be received, acquired, possessed and used only by physicians, clinical laboratories or hospitals and only for *in vitro* clinical or laboratory tests not involving internal or external administration of the material, or the radiation therefrom, to human beings or animals. Its receipt, acquisition, possession, use and transfer are subject to the regulations and a general license of the U.S. Nuclear Regulatory Commission or a State with which the Commission has entered into an agreement for the exercise of regulatory authority.

MP Biomedicals, LLC

Adherence to the basic rules of radiation safety should provide adequate protection. The user is referred to National Bureau of Standards Handbook No. 92, "Safe Handling of Radioactive Materials", issued March 9, 1964, Superintendent of Documents, U.S. Government Printing Office, Washington, D.C. 20402. A summary follows:

♦ Do not eat, drink, smoke or apply cosmetics where radioactive materials are used. ♦ Do not pipet radioactive solutions by mouth. ♦ Avoid direct contact with all radioactive materials by using protective articles such as lab coats and disposable gloves. ♦ All radiological work should be done in a designated area away from traffic. ♦ Radioactive materials should be stored in their original containers in a designated area. ♦ A record book for logging receipt and disposal of all radioactive materials should be kept. ♦ Laboratory equipment and glassware which are subject to contamination should be segregated to prevent cross-contamination of different radioisotopes. ♦ Any radioactive spills should be taken care of immediately in accordance with established procedures. ♦ All radioactive materials must be disposed of in accordance with the prevailing regulations and guidelines of the agencies

3

holding jurisdiction over the laboratory. ♦ Uncontaminated containers may be discarded in non-radioactive waste providing that labels and labeling are defaced.

G. Equipment and Reagents Required but Not Provided in Kit:

Evacuated glass tube, 5 or 10 mL or Evacuated glass tube containing EDTA, 7 or 10 mL.

Polypropylene or polystyrene tubes (12 x 75 mm), disposable.

Test tube rack.

Any type semi-automatic pipette with disposable tips capable of delivering 100 µL, 200 µL and 1.0 mL.

Vortex mixer.

Centrifuge capable of achieving an RCF of at least 1000 x g.

Gamma counter for measuring [<sup>125</sup>I] and [<sup>57</sup>Co] simultaneously or sequentially with adjustable counter settings capable of separating the two counting spectra.

H. Additional Equipment and Reagents Required for Red Cell Folate Assay:

- Ascorbic acid, MP Biomedicals, Catalog No. 06B259110.
- Micro-Hematocrit Centrifuge.
- Micro-Hematocrit Capillary Tubes/Heparinized.
- Holder for sealing and carrying microcapillary tubes.

IV. SPECIMEN COLLECTION:

Samples should be collected from fasting individuals, since recent food intake may increase the folic acid level appreciably. The laboratory should be advised of possible radioactivity in the sample.

Samples must not be collected in ascorbic acid or in high concentrations of fluoride, since either of these two agents appears to destroy vitamin B<sub>12</sub>.

Do not use hemolyzed specimens for serum or plasma assays.

A. Preparation of specimen for analysis:

1. **Serum or Plasma**  
Collect the blood in a 5 or 10 mL evacuated glass tube. If serum is being collected, allow the blood to clot at room temperature in the closed tube for 30-60 minutes. Use EDTA if plasma is desired for analysis. Centrifuge for 10 minutes and collect the serum or plasma.
2. **Whole Blood Hemolysate (for red cell folate assay)**  
Collect the blood in a 7 mL evacuated glass tube containing EDTA. Determine and record the hematocrit value.  
Add 100 µL well-suspended blood to 2 mL of freshly prepared 0.2% ascorbic acid solution (w/v). This is a 1:21 dilution. Mix by inversion several times; avoid foaming.  
Let the hemolysate stand at +20 to +25°C for 60 to 90 minutes prior to assay. Protect from light during this time.

4

*Sample received via evacuated tube 11/2/84*

*[Handwritten signatures and initials]*

**MP BIOMEDICALS**

**LABORATORY**

*[Handwritten signature]*

- B. Shipping of Specimens:  
Carefully packed serum and plasma should be shipped and received frozen.
- C. Storage of Samples:  
Store samples before analysis at 2-8°C. If storage is expected to exceed 4 hours, the sample should be stored at -20°C or below. Samples are stable for 6-8 weeks at this temperature. Do not store in freezers with an automatic defrost cycle in order to avoid repeated freezing and thawing.

V. PROCEDURE:

If the counter has two or more channels it should be calibrated to count [<sup>125</sup>I] in one channel and [<sup>57</sup>Co] in another channel. If the counter has one channel it should be calibrated so that different counter settings will count one isotope at a time. In the case of the latter, it will be necessary to count the tubes twice, once setting the counter for [<sup>125</sup>I] and obtaining the folate curve and sample data, then setting the counter for [<sup>57</sup>Co] to obtain the vitamin B<sub>12</sub> curve and sample data. All reagents and samples must be brought to room temperature before use. Return to recommended storage temperature after use. Do not use reagents other than those provided in this particular kit. Vitamin B<sub>12</sub> and folates are light sensitive and should be exposed to diffused light only, for as short a period of time as possible. It is advisable to cover the rack of assay tubes during the extraction and binding steps.

In the following protocol it is recommended that the standard level points be run in duplicate. Patient samples must be assayed in duplicate and the preparation of the standard curve and the clinical determinations must be run simultaneously. Control sera should be run at the same time as patient samples.

- A. Reagent Preparation:
- If one bottle of tracer can be used within 30 days, add the contents of one vial of DTT to a bottle of tracer. If the same bottle of tracer is to be used for periods longer than 30 days, follow the directions below.
  - Equal amounts (100 µL each) of tracer and dithiothreitol (DTT) are added per assay tube. Mix tracer and dithiothreitol 1 to 1 and use 200 µL/tube. Alternatively the tracer and DTT may be pipetted separately, 100 µL/tube of each, pipetting the tracer first followed by the dithiothreitol.
- B. Assay Procedure:
- Number 16 tubes for the standards. Beginning with 17, number two tubes for each clinical sample.
  - Add standards and clinical samples according to the outline which follows.
  - Add 200 µL Working Tracer/ DTT Solution (Reagent B1) to all tubes including the Total Count tubes (1 and 2). Vortex.
  - Incubate for 15 minutes at room temperature. (18 - 25°C)
  - Add 100 µL Extracting Reagent to tubes 3-16 and all sample tubes. Vortex.
  - Incubate for 10 minutes at room temperature. (18 - 25°C)
  - Thoroughly mix the bottle of SimulTRAC-SNB Blank Reagent. Add 1000 µL blank reagent to tubes 3 and 4.
  - Thoroughly mix the bottle of SimulTRAC-SNB Binder. **Caution: Must be shaken vigorously.** Add 1000 µL binder to tubes 5-16 and all sample tubes. Vortex.
  - Incubate tubes 3-16 and all sample tubes at room temperature (18 - 25°C) for 60 minutes from the time of the last addition of the binder. Cover the rack of tubes with aluminum foil to exclude light, or keep in

a dark location.

- Centrifuge at a minimum of 1000 x g for 10 minutes, preferably in the cold.
- Gently decant and discard each supernatant. Remove the last drop by touching the tube to a paper towel or absorbent paper.
- Count the radioactivity in the pellets and in tubes 1 and 2 in sequence for one minute with a gamma counter. The total count per minute for tubes 1 and 2 for [<sup>57</sup>Co] should be between 10,000 and 25,000 and for [<sup>125</sup>I] between 15,000 and 35,000, depending on the instrument and age of the tracer. A shorter counting time may be used provided the counts in tubes 1 and 2 are at least 10,000 (total count).

SimulTRAC-SNB RADIOASSAY

Tube	Standard or Sample (µL)	Working Tracer Solution (µL)	Incubate	Extracting Reagent (µL)	Incubate	Blank Reagent (µL)	Binder (µL)	Incubate	Centrifuge
1, 2	---	200	Vortex.	---	Vortex.	---	---	Vortex.	Centrifuge
3, 4	200A	200	Incubate	100	Incubate	1000	---	Incubate	all tubes
5, 6	200A	200	all tubes	100	all tubes	---	1000	all tubes	(except 1
7, 8	200B	200	at room	100	at room	---	1000	at room	and 2) at
9, 10	200C	200	temper-	100	temper-	---	1000	temper-	1000 x g for
11, 12	200D	200	ature (18	100	ature (18	---	1000	ature (18	10 min.
13, 14	200E	200	- 25°C)	100	- 25°C)	---	1000	- 25°C)	
15, 16	200F	200	for 15	100	for 10	---	1000	for 60	
min.					min.			min.	
Patient Samples	200	200		100		---	1000		

After centrifugation, decant the supernatants and count the radioactivity in the pellets.

Calculate. Draw standard curve and determine patient assay values.

VI. CALCULATION OF RESULTS:

A. Serum or Plasma Vitamin B<sub>12</sub> and Folate:

The Vitamin B<sub>12</sub> curve and values are calculated from the data obtained by counting [<sup>57</sup>Co]; the folate curve and values are calculated from the [<sup>125</sup>I] counts.

1. Average the counts found in tubes 3 and 4, the "Blank" tubes. Subtract the "Blank" from all other tube counts to obtain the corrected counts. Use only the corrected counts in the calculations. **NOTE:** The unit of time must be constant for all tubes counted.
2. Average the corrected counts for tubes 1 and 2 to give the corrected Total Count per assay.
3. Divide the average of the corrected counts for tubes 5 and 6 by the corrected Total Count to give the Trace Binding, B<sub>0</sub>. This value should be greater than 35%.
4. Divide the corrected counts for each tube by the average corrected counts for tubes 5 and 6 to give the % of Trace Binding for each tube. A Standard Curve may be plotted as follows:  
Using logit-log paper, plot % of Trace Binding as the ordinate versus pg/mL vitamin B<sub>12</sub> or ng/mL folate standard on the log scale. Typical counting data and calculated % of Trace Binding are given in Table 1; the standard curves plotted for these data are shown in Figure 1. In practice, it may be advisable to plot each standard curve on a separate sheet to avoid possible confusion.
6. The concentration of vitamin B<sub>12</sub> or folate in serum or plasma is determined by interpolation from the standard curve of % of Trace Binding versus either pg/mL vitamin B<sub>12</sub> or ng/mL folate (Figure 1).
7. Example:
  - a. The following example is for a B<sub>12</sub> sample. The same calculation procedure is used for folate:

$$\text{Blank} = \frac{\text{counts tube 3} + \text{counts tube 4}}{2} = \frac{750 + 728}{2} = 739$$

$$\text{Corrected Total Count} = \frac{\text{Counts (tube 1-Blank)} + (\text{tube 2-Blank})}{2} = \frac{(23319 - 739) + (23716 - 739)}{2} = 22778$$

$$B_0 = \text{Trace Binding} = \frac{\text{av. corr. counts for tubes 5 \& 6}}{\text{corrected Total count}} \times 100 = \frac{10632}{22778} \times 100 = 46.7\%$$

% of Trace Binding

$$\text{calculation for Tube 7} = \frac{(\text{counts tube 7} - \text{Blank})}{\text{av. corr. counts for tubes 5 \& 6}} \times 100 = \frac{(9914 - 739)}{10632} \times 100 = 86.2\%$$

Sample Calculation for a Patient Specimen:

$$\text{Count (found)} = 7324$$

$$\text{Blank} = 739$$

$$\% \text{ of Trace Binding} = \frac{7324 - 739}{10632} \times 100 = 61.9\%$$

The logit-log Standard Curve (Figure 1) shows that 61.9% corresponds to a B<sub>12</sub> concentration of 470 pg/mL.

The average of this value and the value for the duplicate determination is reported as the B<sub>12</sub> concentration in pg/mL for the patient sample.

B. Red Cell Folate:

1. Calculate the % of Trace Binding for the hemolysate (Step A-4)
2. Obtain the concentration of folate by interpolation from the Standard Curve (Step A-6).
3. Multiply the folate concentration by 21. (A 1:21 dilution of whole blood was made in preparing the specimen.) This gives the folate concentration in ng/mL of whole blood.
4. Divide the folate concentration of whole blood by the hematocrit expressed as a decimal. This gives the folate concentration in ng/mL of packed red cells.

$$\text{ng/mL of packed red cells} = \frac{(\text{ng/mL of hemolysate}) \times 21}{\% \text{ hematocrit}/100}$$

Sample Calculation for Red Cell Folate:

$$\text{Count (found)} = 11271$$

$$\text{Blank} = 1321$$

$$\% \text{ of Trace Binding} = 56.3$$

The logit-log Standard Curve (Figure 1) shows that 56.3% corresponds to a folate concentration of 4.6 ng/mL.

$$\text{Patient hematocrit} = 42\%$$

$$\text{ng/mL of packed red cells} = \frac{4.6 \times 21}{0.42} = 230 \text{ ng/mL}$$

The average of this value and the value for the duplicate determination is reported as the red cell folate concentration in ng/mL for the patient sample.

For serum vitamin B<sub>12</sub> ("true" cobalamin) expected values, two populations were used. The first population consisted of 38 vitamin B<sub>12</sub> deficient subjects with confirmed diagnoses including pernicious anemia, gastric or intestinal damage, or disease states associated with B<sub>12</sub> deficiency. The estimated 99 percentile of this population defined the upper limit of the deficient (low) range.

The second population consisted of 121 healthy male and female volunteers ranging in age from 19 to 60, with no unusual dietary habits. These volunteers were judged to be hematologically normal by standard laboratory criteria. The central 95 percent of this population defined the expected range (normal). The indeterminate range is defined as the range of values between deficient and normal populations.

Interpretation	Vitamin B <sub>12</sub> pg/mL	(pmol/L)
Low	<120	<88
Indeterminate	120 - 160	88 - 118
Normal	160 - 970	118 - 716
High	>970	>716

In another study specimens were compared by the SimulTRAC-SNB assay and the *Euglena gracilis* microbiological assay. This is shown in Figure 2.

Serum vitamin B<sub>12</sub> levels well above 1000 pg/mL are suggestive of either liver disease or a myeloproliferative disorder such as polycythemia vera, myeloid metaplasia or chronic granulocytic leukemia. Levels above 4000 pg/mL are unusual in liver disease, but are common in myeloproliferative disorders<sup>6,13,14</sup>.

IX. SPECIFIC PERFORMANCE CHARACTERISTICS:

A. Accuracy:

1. A comparison of the results obtained for vitamin B<sub>12</sub> with this kit and those obtained using the SimulTRAC-S Radioassay kit gave these regression data:

$$(x = \text{SimulTRAC-S}, y = \text{SimulTRAC-SNB})$$

Number of samples = 160      Correlation coefficient = 0.98  
Slope = 0.93      Y-Intercept = -26 pg/mL

2. Comparison of results with those obtained by microbiological assay with *Euglena gracilis*<sup>15</sup> gave the following calculated regression.

$$(x = \text{microbiological assay}, y = \text{SimulTRAC-SNB}):$$

Number of samples = 132      Correlation coefficient = 0.96  
Slope = 1.1      Y-Intercept = -32 pg/mL

3. A comparison of the serum folate results obtained in this kit with those obtained using the SimulTRAC-S kit gave the following results:

Number of samples = 160      Correlation coefficient = 0.97  
Slope = 1.03      Y-Intercept = -0.01 ng/mL

B. Precision:

1. Intra-assay variation:

a. Vitamin B<sub>12</sub>

Specimen	N	Mean pg/mL	S.D.	% C.V.
Control 1	20	318	19.5	6.1
Control 2	22	444	14.1	3.2
Control 3	20	593	26.3	4.4
Control 4	20	1203	69	5.7
Serum Pool	20	169	19	11.2

b. Folate

Specimen	N	Mean ng/mL	S.D.	% C.V.
Control 1	20	1.87	0.16	8.6
Control 2	22	2.16	0.15	6.9
Control 3	20	4.89	0.25	5.1
Control 4	20	12.1	0.54	4.5
Serum Pool	20	6.10	0.25	4.1

2. Inter-assay variation:

a. Vitamin B<sub>12</sub>

Specimen	N	Mean pg/mL	S.D.	% C.V.
Control 1	35	313	25.6	8.2
Control 2	44	432	27.9	6.4
Control 3	52	555	38	6.8
Control 4	52	1106	47	4.2
Serum Pool	30	151	18.6	12.3

b. Folate

Specimen	N	Mean ng/mL	S.D.	% C.V.
Control 1	35	1.62	0.19	11.7
Control 2	44	1.96	0.16	8.2
Control 3	52	4.50	0.34	7.5
Control 4	52	10.4	1.0	9.6
Serum Pool	30	6.60	0.47	7.1

Handwritten signature and notes at the bottom of page 12, including a circular stamp and the text "W. A. R. 2 8 0 120".

C. Recovery:

One serum sample was spiked with cyanocobalamin and one with MTFA. The MTFA was calibrated by UV spectrophotometry and spiked into the sample. Recoveries were calculated as:

$$\frac{\text{found} - \text{endogenous}}{\text{spike}} \times 100$$

Vitamin B <sub>12</sub>				MTFA			
Spike	Found	Expected	% Recovery	Spike	Found	Expected	% Recovery
0	340	---	---	0	0.9	---	---
200	541	540	100	2	3.03	2.9	107
400	763	740	106	4	4.95	4.9	101
1000	1297	1340	96	10	10.66	10.9	98

D. Linearity

A serum sample was diluted with standard A to test for linearity. Linearity was calculated as:

$$\frac{\text{found}}{\text{expected}} \times 100$$

Vitamin B <sub>12</sub>				Folate			
Dilution	Found	Expected	%	Dilution	Found	Expected	%
0	398	---	---	0	19.48	---	---
1.33	273	299	91	1.33	14.40	14.6	99
2	195	199	98	2	10.17	9.74	104
4	100	100	101	4	5.49	4.87	113
8	50.0	50	100	8	2.68	2.44	110

E. Specificity

5-Methyltetrahydrofolic acid and pteroylglutamic acid have equal affinity for the binder in this assay.

The porcine Intrinsic Factor used in this kit has been purified by affinity chromatography and contains less than 4% R protein. The purity has been established according to criteria proposed by the National Committee for Clinical Laboratory Standards<sup>18</sup>.

The intrinsic factor in the MP Biomedicals SimuTRAC-SNB Radioassay Kit measures 10,000 pg/mL of cobinamide as less than 75 pg/mL when read against a standard curve constructed from cobalamin standards. In addition, vitamin B<sub>12</sub> binding is more than 95% inhibited by specific anti-IF blocking antibody.

Additional studies have shown that when using the MP Biomedicals SimuTRAC-SNB Kit, cobinamide (a vitamin B<sub>12</sub> analogue) is less than 0.01% cross-reactive.

F. Sensitivity:

Sensitivity as defined by the concentration at 90% trace binding is 75 pg/mL for B<sub>12</sub> and 0.6 ng/mL for folate.

Handwritten signature and stamp at the bottom of page 14. The stamp is circular and contains the text "VITAMIN B12" and "COBALAMIN".

## SimulTRAC-SNB Radioassay Kit Vitamin B<sub>12</sub> [<sup>57</sup>Co]/Folsäure [<sup>125</sup>J]

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Figure 1 - SimuTRAC-SNB  
Typical Standard Curve with Logit-Log Transformation



Figure 2

*[Handwritten signatures and notes]*

*[Circular stamp]*

*[Handwritten text: "Dr. Brewer", "10/20/80", "D.L. 80A"]*

**BIBLIOGRAPHY**

1. Herbert, V., Introduction to the "Nutritional Anemias". Sem. Hemat., 7, 2-5, 1970.
2. Herbert, V., Folic Acid and Vitamin B<sub>12</sub>. "Modern Nutrition in Health and Disease", 5th ed., edited by R.S. Goodhart and M.E. Shils, Lea & Febiger, Philadelphia, pp. 221-224, 1973.
3. Chanarin, I., The assay and concentration of folate in blood and other tissues. "The Megaloblastic Anemias", Blackwell Scientific Publications, Oxford, pp. 306-336, 1969.
4. Herbert, V., Drugs effective in megaloblastic anemias. "The Pharmacological Basis of Therapeutics", 4th ed., edited by L.S. Goodman and A. Gilman, MacMillan Co., New York, pp. 1414-1444, 1970.
5. Tishman, G. and Herbert, V., B<sub>12</sub> dependence of cell uptake of serum folate: an explanation for high serum folate and cell folate depletion in B<sub>12</sub> deficiency. Blood, 41, 465-469, 1973.
6. Herbert, V., B<sub>12</sub> and folate deficiency. "Nuclear Medicine", edited by B. Rothfield, Lippincott, Philadelphia, pp. 69-84, 1974.
7. Chanarin, I., The assay and concentration of Vitamin B<sub>12</sub> in serum, liver and cerebrospinal fluid. "The Megaloblastic Anaemias", Blackwell Scientific Publications, Oxford, pp. 192-229, 1969.
8. Sullivan, L.W., The megaloblastic anemias. "Hematology: Principles and Practice", edited by E.C. Mengel, E. Frei III and R. Nachman, Yearbook Medical Publishers, Chicago, pp. 95-131, 1972.
9. Kohouse, J.F., Kondo, H., Allen, N.C., Podell, E. and Allen, R.H., Cobalamin analogues are present in human plasma and can mask cobalamin deficiency because current radioisotope dilution assays are not specific for true cobalamin. New Eng. J. Med., 299, 785-792, 1978.
10. Cooper, B.A. and Whitehead, V.M., Evidence that some patients with pernicious anemia are not recognized by radiodilution assay for cobalamin in serum. New Eng. J. Med., 299, 816-818, 1978.
11. Zucker, R.M., Podell, E.R., and Allen R.H., Multiple Problems with Current No Boil assays for Serum Cobalamin, Ligand Quarterly, 4, 54-58, 1981.
12. Higgins, T., and Wu, A., Differences in Vitamin B<sub>12</sub> Results as Measured with Boil and No Boil Kits, Clin. Chem., 29, 587-588, 1983.
13. Herbert, V., Diagnostic and prognostic values of measurement of serum Vitamin B<sub>12</sub>-binding proteins. Blood, 32, 305-312, 1968.
14. Gilbert, H.S., Krauss, S., Pasternack, B., Herbert, V. and Wasserman, L.R., Serum Vitamin B<sub>12</sub> and unsaturated Vitamin B<sub>12</sub>-binding capacity in myeloproliferative disease: value in differential diagnosis and as indicators of disease activity. Ann. Intern. Med., 71, 719-729, 1969.
15. Anderson, B., Investigations into the Euglena method for the assay of Vitamin B<sub>12</sub> in serum. J. Clin. Pathol., 17, 14-26, 1964.
16. Guidelines for Evaluating a B<sub>12</sub> (Cobalamin) Assay. National Committee for Clinical Laboratory Standards, March 1980.

**Manufactured by:**  
**Hergestellt von:**  
**Fabriqué par:**  
**Prodotto da:**  
**Fabricado por:**

**MP Biomedicals, LLC**  
**Diagnostics Division**  
**29525 Fountain Parkway**  
**Solon, Ohio 44139**  
**USA**

**Customer Service:**  
**FAX:**  
**Technical Service:**

**(800) 854-0530**  
**(440) 337-1180**  
**(800) 854-0530**  
[www.mpbio.com](http://www.mpbio.com)

**European**  
**Authorized**  
**Representative:**

**MP Biomedicals Germany GmbH**  
**Thüringer Straße 15**  
**37269 Eschwege**  
**Tel: +49 5651 9210**  
**Fax: +49 5651 921 181**

CE

06B58581-R16  
 Q15-042, A15-005  
 Q16-009, A16-001  
 04/16

*[Handwritten signatures and stamps]*

*[Circular stamp: "MP BIOMEDICALS" with illegible text]*

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Всички права запазени © 2018 ИЗПЪЛНИТЕЛНА АГЕНЦИЯ ПО ЛЕКАРСТВАТА